



**DISPOSABLE &
ACCESSORIES**



Made in Tuscany with care



SALIVA EJECTORS

p.5



DISPOSABLE TOWELS

p.13



ACCESSORIES

p.21



Over 100 years of history

A family story combining expertise and distinguishing Italian quality searching for technologically advanced, innovative, customized and flexible solutions that meet the highest standards demanded by the international markets.

Since 1987 we have developed and manufactured solutions for the dentistry, medical and industrial markets.

A stable and well-structured company being part of the Nostrum Group that has manufactured plastic and dentistry items since 1975.

We merge the special care of Italian handcraft with constant research into technological and innovative solutions with the sole aim of meeting the demands of the market.



1923

The Nostrum Group was founded in 1923 as manufacturer of automotive gas generators with a company branch skilled in moulding of plastic components.

1975 - 2025 50 YEARS PRODUCING

50 years of history since the Grotti family took over the Nostrum group as the main shareholder. It started with galvanic and anodizing process while strengthening the plastic moulding.

In 1987 Lucia Puccetti and Luciano Grotti established Dental Market, a Company specialized in the production of accessories for dental surgery with 3 primary components made of plastic, stainless steel and anodized aluminium.


2018

Grotti's daughters pursue the group path taking up a new challenge as female entrepreneur by renewing the Company and the Brand PLUSS production, with passion, to chase excellence, quality and innovation with a new and distinct vocation to international markets.

PLUSS production mission is to be a dental and medical atelier, focused on Made in Italy products. Our aim is to satisfy any production demand being very versatile both as to quantity and design.

Our technical department can fully develop any customized projects in terms of quantities and design thanks to our main strengths:

- Design and development production of metal - stainless steel and aluminium items
- Moulding and extrusion of plastic raw materials
- Production of disposable tissue articles



Quality Made in Italy

PLUSS Made in Italy means quality, high specialisation, distinction and style. We are able to add an handmade touch to productions.




We believe in a healthier and more sustainable economy and our values are:

- Circular economy
- Efficient choice and use of raw materials
- Development of with long shelf life quality products
- Patented no glue policy
- Research of compostable and recycled raw materials
- Recycling practices

We believe to raise awareness among final consumers because every action have a good impact for the Planet.

Our production model

A production method springing from decades of experience allows us to have always any items in stock ready for dispatch, thus avoiding you needless stockpiling and, at the same time, ensuring the best value for money while still giving the chance of receiving even small lots of product.



Our abiding goal

is merging our long-term experience of activity with production flexibility, keeping up with the pace of the market's innovations.

We ceaselessly innovate to find new solutions to your demands looking for the highest production reliability by optimising our processes.



SALIVA EJECTORS

Made in Tuscany with care





Patented Quality

The continuous research and improvements led to achieving three patents PLUSS Made in Italy saliva ejectors.

1

First Patent



The inner wire is fastened to the tube without glue, as other manufacturers do. It is fused into the plastic during the production process.

The result is a unique and entirely non-toxic product. This method ensures that the wire not come out when bending the ejector



2

Second Patent



More is Less. The 100-piece bag is fastened through a safety band entirely wrapping around it. An extra advantage we grant to dentists over our competitors. By fixing the band, we assure our customers of carrying out all the quality controls during production. It is also a transport cost saving for the retailers since the carton containing the wrapped bags is smaller than the standard packaging. A 120×80 pallet can hold nearly 50% more goods than the standard one without damaging the lowest cartons and with a clear saving in transport costs.

3

Third Patent



While manufacturing, our saliva ejectors undergo an individual 100-piece bag decontamination process, thus ensuring cleaner products than our competitors. Our saliva ejectors are the only ones with these characteristics on the international market.

Non-Toxic System

Our patented system eliminates any traces of glue or dangerous substances from the saliva ejectors through high-performance and safe production and clamping system.




The only glue free wire

Non toxic fused in wire
Our metal wire is fused into the ejectors during the plastic extrusion process, thus ensuring the highest tube flexibility without risking it unloking.

3D secure cap lock

OUR NON TOXIC IMMUABLE TIP



The tip is mechanically attached to the ejector trough the perfect interlocking of hooks integrated into the tip mould. The ejector does not allow the tip to come off, blocking its position persistently.


LESS IS MORE

Patented

50% more goods than the standard one without damaging the lowest cartons and with a clear saving in transport costs.

The box is considerably smaller compared to the competitor's box. This allowed a clear saving in transport costs.







**DECONTAMINATED
PRODUCT**



NO GLUE POLICY



LONG SALIVA EJECTORS

Lenght 145 mm
Pcs100 - Box 5000



Code	Color
3410/01/1/008	Transparent
3410/12/1/008	Opaque Light Blue
3410/02/1/008	Opaque White
3410/03/1/008	Opaque Green
3410/08/1/008	Orange
3410/09/1/008	Light Blue
3410/04/1/008	Blue
3410/11/1/008	Fucsia
3410/06/1/008	Yellow
3410/13/1/008	Lilac
3410/07/1/008	Red
3410/10/1/008	Pink
3410/05/1/008	Green
3410/14/1/008	Fresh Green




SHORT SALIVA EJECTORS

Lenght 130 mm
Pack 100 - Box 5000



Code	Color
3400/01/1/008	Transparent
3400/12/1/008	Opaque Light Blue
3400/02/1/008	Opaque White
3400/03/1/008	Opaque Green
3400/08/1/008	Orange
3400/04/1/008	Blue
3400/11/1/008	Fuchsia
3400/06/1/008	Yellow
3400/13/1/008	Lilac
3400/07/1/008	Red
3400/10/1/008	Pink
3400/05/1/008	Green
3400/14/1/008	Fresh Green



DISPOSABLE TOWELS

Made in Tuscany with care



Disposable Towels

Made in Tuscany with care.



Our Disposable Towels. A Made in Italy value-added

Made of 1 or 2 pure cellulose an recycled plies, they ensure a perfect absorption of liquids leaking out during treatments, as well as ideal impermeability, thanks to the polyester film effectively laminated without adhesives.

Manufactured in our facility in Lucca, the beating heart of Italian paper production, from top-quality and entirely non-toxic raw materials.

- Dimensions
330 x 450 mm
- Class
Class I medical device
- Package
500-loose dental towels box.
50-dental towels package.
- Colours
15 + different colours available.

- White
- Yellow
- Orange
- Peach Pink
- Pink
- Fuchsia
- Magenta
- Lilac
- Dark blue
- Light blue
- Lime
- Emerald green
- Green
- Black
- Sand

Made in Tuscany with care.



Made in Tuscany with care.



NO GLUE AND NO TOXIC POLICY

In our plant in Lucca, which is worldwide known as one of the leading paper districts, we produce a wide range of tissue items from top quality raw materials in a full assortment of colours and sizes.

- Everything is Made in Italy in our facilities in Tuscany.
- Stock ready to order.
- Small quantity orders and customisable



PLUSS Disposable towels completely non-toxic as they do not contain any type of glue. Made in Italy

- 100% recycled tissue (light blue, green and sand)
- Eco Label Certification for sand, white and light blue colours
- High quality of raw materials

Made in Tuscany with care.



Made in Tuscany with care.



DISPOSABLE TOWELS 3 PLIES**(2 plies paper+1 ply pe)**

Lenght 330x450 mm

Pcs/Box 500

**Code Color**

TD0850AZZURRO-08	Blue
TD0850VERDE-08	Green
TD0850BIANCO-08	White
TD0850GIALLO-08	Yellow
TD0850LILLA-08	Lilac
TD0850PESCA-08	Peach Pink
TD0850ROSA-08	Pink
TD0850NERO-08	Black
TD0850SABBIA-08	Sand
TD0850BLUNOTTE-08	Dark Blue
TD0850LIME-08	Lime

**DISPOSABLE TOWELS 2 PLIES****(1 plies paper+1 ply pe)**

Lenght 330x450 mm

Pcs/Box 500

**Code Color**


TD0851AZZURRO-08	Blue
TD0851VERDE-08	Green
TD0851BIANCO-08	White
TD0851ARANCIONE-08	Orange
TD0851FUCCIA-08	Fuchsia
TD0851MAGENTA-08	Magenta
TD0851SMERALDO-08	Emerald
TD0851BLUNOTTE-08	Dark Blue



SMALL STEPS, FOR A BIG IMPACT



 100% compostable

 Eco-friendly

 Absorbent and Durable

 **PLUSS**[®]

COMPOSTABLE TOWELS 3 PLIES**(2 plies paper+1 ply compostable film)**

Lenght 330x450 mm

Pcs/Box 500



2 PLY PAPER 100% RECYCLED WITHOUT CHEMICAL ADDITIVES SAND COLOR







1 PLY FILM: BIOPLASTIC, BIODEGRADABLE, COMPOSTABLE

Code**Color**

TD0850ECO

Sand

**OUR PLUSS:**

-  **RECYCLED TISSUE:** Sand color, 100% free of Chemical Additives with Eco- Label Certification
-  **COMPOSTABLE FILM:** Compostable product conforming to UNI EN ISO 13432, Certificate "Ok Compost"
-  It contains a part of monomers of bio and vegetable origin
-  No glue
-  No toxic policy
-  Produced in Italy

ACCESSORIES

Made in Tuscany with care





STERILIZABLE RUBBER BIBS HOLDERS - mixed color

TD0811 - Pcs/Pack 4



STERILIZABLE BIBS HOLDERS - color Silver

TD0805 - Pcs/Pack 1

DISPOSABLE DAPPEN BIG

50 pcs.



Code	Color	
TD0650	Red	
TD0651	Blue	
TD0652	Green	
TD0653	Yellow	
TD0654	White	
TD0655	Mix	200 pcs

DISPOSABLE PROTECTIVE SHIELD

2 pcs.



TD0700 Shield with elastic

COTTON PELLETS

6 pcs./ Box



TD0642	small	ø 4mm
TD0643	medium	ø 5.5mm
TD0644	large	ø 7mm

CALCINABLE PINS

60 pcs. - Length 23mm



TD7055	Red
TD7056	Green

CALCINABLE PINS

100 pcs.



TD7050	Yellow	n°0/0.6 mm
TD7051	Blue	n°1/0.8 mm
TD7052	Transparent	n°2/1 mm
TD7053	Transparent	n°1/0.8 mm


INSTANTANEOUSLY DISPOSABLE ICE

25 pcs.



TD0900	disposable tnt instant ice
TD0901	disposable polythene instant ice

STERILIZATION FLAT REELS



TD9500	55x200 mt.
TD9501	75x200 mt.
TD9502	100x200 mt.
TD9503	150x200 mt.
TD9504	200x200 mt.
TD9505	250x200 mt.

STERILIZATION SELF SEAL POUCHES

200 pcs.



TD9520	90x230 mm	lung. tot. 260mm
TD9521	140x250 mm	lung.tot. 295mm

DISPOSABLE WASTE RECEIVER

1000 pcs.



TD3065 - White disposable waste receiver

ALLUMINIUM SUPPORT FOR WASTE RECEIVER



TD3066ARG		Silver
TD3066BLU		Blue
TD3066ORO		Gold
TD3066NERO		Black

PLASTIC MIXING PAD



TD7090	50x50 mm	6 pcs.
TD7091	80x80 mm	6 pcs.
TD7092	140x80 mm	3 pcs.

PAPER PLASTIFIED MIXING PAD



TD7095	50x50 mm - 6 pcs.	6 pcs.
TD7096	100x100 mm - 6 pcs.	6 pcs.
TD7097	120x180 mm - 3 pcs.	3 pcs.

POLYETHYLENE AND PAPER DISPOSABLE BIBS

50x60 cm - 80 tears



TD0855ARANCIONE	Orange
TD0855BIANCO	White
TD0855BLU	Blue
TD0855BLUMED	Dark Blue
TD0855LILLA	Lilac
TD0855VERDE	Green

POLYSTYRENE DISPOSABLE LARGE TRAY

280x183x17 mm pcs/box 400



TD3060 smooth surface

TD3060Z ribbed surface

POLYSTYRENE DISPOSABLE MINI TRAY

140x183x17mm pcs/box 400



TD3061 smooth surface

TD3061Z ribbed surface

POLYSTYRENE DISPOSABLE MAXI BASIN

200x150 mm pcs/box 400



TD3062 smooth surface

TD3062Z ribbed surface

POLYSTYRENE DISPOSABLE MINI BASIN

200x100 mm pcs/box 400



TD3063 smooth surface

TD3063Z ribbed surface

DISPOSABLE MIXING PALLET AND SPATULA



TD0338	Spatula - pcs/box 100
TD3064	Mixing pallet - pcs/box 1000

INTERDENTAL PLASTIC WEDGES

kit 100 wedges



TD7060	Yellow	XS
TD7061	Green	S
TD7062	Orange	M
TD7063	Lilac	L
TD7064	Mix	All

MIXING BOWLS FOR PLASTER AND ALGINAT

GREEN



TD0321	Ø 100mm
TD0322	Ø 120mm
TD0323	Ø 140mm
TD0324	Ø 160mm

WHITE



TD0331	Ø 100mm
TD0332	Ø 120mm
TD0333	Ø 140mm
TD0334	Ø 160mm

FLEXIBLE SPATULA



TD0335	flexible spatula	Light Blue
TD0336	hard spatula	Blue
TD0337	flexible anatomic spatula	White

COLOURED FLEXIBLE SPATULA



TD0337MIX	flexible spatula	6 Colors
-----------	------------------	----------

AUTOCCLAVABLE PLASTIC COTTON PELLETS DISPENSER

50x50x53 mm



TD0645	Small
TD0646	Medium
TD0647	Large

METAL CUPS DISPENSER

85x85x350 mm



DS0625	Stainless Steel
DS0626/02	White
DS0626/04	Blue
DS0626/05	Green
DS0626/06	Yellow
DS0626/07	Red
DS0626/08	Orange
DS0626/09	Light Blue
DS0626/13	Lilac
DS0626/14	Lime
DS0626/17	Dark Blue
DS0626/18	Emerald
DS0626/19	Black
DS0626/22	Petrol

METAL GLOVES DISPENSER

230x115x50mm



DS0635	Stainless Steel
DS0635/02	White
DS0635/19	Black
DS0635/23	Gold

METAL DISPOSABLE BIBS DISPENSER

275x125x130 mm



DS0673	Stainless Steel
DS0673/02	White
DS0673/19	Black
DS0673/23	Gold

METAL ICE DISPENSER

180x145x140 mm




DS0690	Stainless Steel
DS0690/02	White
DS0690/19	Black
DS0690/23	Gold

STAINLESS STEEL TOWELS DISPENSER

Front opening with cover






DS0673		Stainless Steel
DS0673/02		White
DS0673/19		Black
DS0673/23		Gold

STAINLESS STEEL MASK DISPENSER

Front opening with cover



DS0671		Stainless Steel
DS0671/02		White
DS0671/19		Black
DS0671/23		Gold

COTTON GAUZE DISPENSER

52x52x55 mm



DS0679

METAL COTTON PELLETS DISPENSER



DS0641

DS0647

DS0641	40x40x50 mm
DS0647	52x52x55 mm

ALLUMINIUM WASTE RECEIVER

Ø 65 h 70 mm



TD3240ARG		Silver
TD3240BIA		White
TD3240BLU		Blue
TD3240ORO		Gold
TD3240ROS		Red
TD3240VER		Green
TD3260		Silver mini Ø 50 H 40 mm

ALLUMINIUM COTTON DISPENSER

Ø 50 h 70 mm



TD3250ARG		Silver
TD3250BIA		White
TD3250BLU		Blue
TD3250ORO		Gold
TD3250ROS		Red
TD3250VER		Green
TD3265		Silver mini Ø 50 H 40 mm

OMNITRAY FOR STERILIZATION



TD3100 1lt White 260x110x80mm

TD3110 3lt White 350x260x140mm

DOUBLE PLASTIC DISPENSER FOR GAUZES

110x55x80 mm



TD0689

for gauzes 5x5 cm

PLEXIGLASS MASKS DISPENSER

200x120x110 mm



TD0671

transparent with cover

PLASTIC TOWELS DISPENSER

180x130x170 mm




TD672BIA

front opening

PLASTIC BIBS DISPENSER

280x136x125 mm



TD0674BIA

dispenser for bibs in rolls



 **PLUSS**[®]

PLASTIC ORTHODONTIC BOX

65x76x40 mm - 10 pcs/pack



TD7021	Mix
TD7021ARANCIONE	Orange
TD7021AZZURRO	Light blue
TD7021BIANCO	White
TD7021BLU	Blue
TD7021FUCSIA	Fuchsia
TD7021GIALLO	Yellow
TD7021GIALLOFLUO	Fluo yellow
TD7021GRIGIO	Grey
TD7021NERO	Black
TD7021ORO	Gold
TD7021ROSA	Pink
TD7021ROSSO	Red
TD7021SMERALDO	Emerald
TD7021VERDE	Green
TD7021VERDECHIARO	Fresh green
TD7021VIOLA	Purple



BOX FOR BITE

78x76x30 mm - 10 pcs/pack



TD7022ARANCIONE	Orange
TD7022AZZURRO	Light blue
TD7022BIANCO	White
TD7022BLU	Blue
TD7022FUCSIA	Fuchsia
TD7022GIALLO	Yellow
TD7022LILLA	Lilac
TD7022LIME	Lime
TD7022ROSSO	Red
TD7022VERDE	Green
TD7022BLUMETAL	Glittered blue

BOX FOR PROSTHESIS WITH FILTER

Ø90 h.82 mm - 10 pcs/pack



TD7025GIALLO	Yellow
TD7025VERDE	Green
TD7025AZZURRO	Light blue
TD7025BIANCO	White
TD7025BLU	Blue
TD7025FUCSIA	Fuchsia

PLASTIC ORTHODONTIC BOX

55x75x40 mm - 10 pcs/pack



TD7020	Mix
--------	-----

PLEXIGLASS SALIVA EJECTOR DISPENSER

160x140x115mm



TD0670 Transparent

PLASTIC COTTON ROLLS DISPENSER

TD0677 1 compartment H 100 Ø100mm

TD0678 4 compartment H 50 Ø124m

PLASTIC COTTON ROLLS DISPENSER

160x60x150 mm



TD0676BIA White

TD0676BLU Blue

PLASTIC COTTON ROLLS DISPENSER

95x50x100 mm



TD0675BIA White

TD0675BLU Blue

TD0675GIA Yellow

TD0675GRI Grey

TD0675ROS Red

DESIGNED FOR YOUR THOUGHTS

A space for you to reflect, to create, to plan.
Dentistry is not just about precision, but also about
innovation and vision.

This image shows a single sheet of white paper with horizontal blue lines, similar to standard notebook paper. The lines are evenly spaced and run across the width of the page. There is no handwriting or other markings on the paper.



EXCLUSIVE ITALIAN SOLUTIONS FOR PROFESSIONALS

**Lucca - Italy
+ 39 0583 1900 385
info@dentalmarket.biz
www.plussproduction.it**

